

# Incident Report: Suspicious File Detected in MicroDicom Viewer

## 1. Source of the File

The installer files were downloaded from the official MicroDicom website:  
<https://www.microdicom.com>.

## 2. Filename & Hashes

File	SHA256 Hash	Description	File Size
MicroDicom-2025.1-x64.exe	15ee7e246846a9c46aac53a855f4e2bc123de8c50cfa39dbffd9e3fcdac9b811	Installer from April 2025	12 MB
MicroDicom-2025.3-x64.exe	783af5c29ed89ae46bc62b5775729d46622e3ebcdba628eb85e3558074b22513	Current official installer	14 MB
~~VisusStart er.bin	058370b70c65c1764960ce0b3eaddcd789b473e90097d247b97dec63959cd4f6	Suspicious temporary binary	1.3 MB

## 3. Date & Time of Initial Detection

The suspicious file was first detected on June 5, 2025, approximately at 00:30.

## 4. Tools Used for Analysis

- Malwarebytes (up to date), performed a full system scan. The suspicious file was detected early during an 8-hour scan.
- Windows Defender (up to date), performed a full system scan (exact scan time not available).

- Additional scanning services used: VirusTotal (virustotal.com), FileScan.io, Dr.Web Virus Monitor (vms.drweb.com), IPQualityScore (ipqualityscore.com).

## 5. Findings & Behavior

### Malwarebytes:

- Detection name: Malware.Sandbox.50
- Detection time: June 5, 2025, 00:30 (approx)
- File path: C:\Users\49157\AppData\Local\Temp\~~VisusStarter.bin
- Type: File

### VirusTotal:

- The file ~~VisusStarter.bin was flagged as malicious by 4 out of 62 vendors, including Avira, GData, Cynet, and WithSecure.

### FileScan.io Highlights (only suspicious/important findings):

- XOR decoding loop detected in entry point (Defense Evasion: Obfuscated Files or Information)
- PE imports APIs related to privilege modification, file attribute modification, code injection, and anti-debugging (Defense Evasion, Discovery, Impact)
- Presence of overlay and suspicious section names
- PE resources exceed 75% of total file size

Other results from VirusTotal, Dr.Web, and IPQualityScore showed no significant findings for the installer files.

## 6. Additional Files Affected

Only the temporary file ~~VisusStarter.bin was found to be suspicious.

The installer files themselves showed fewer signs of malicious behavior, with only FileScan.io flagging the current installer.

## 7. System Environment

Operating System: Windows 10, 64-bit, fully up to date.

Malwarebytes and Windows Defender were both updated to the latest versions before scans.

## 8. Context & Background

This analysis was conducted by **Darius Nasser-Jafar**, a cybersecurity enthusiast preparing to start a formal IT training program in cybersecurity. The initial reason for investigation was to view X-ray images of his dog using MicroDicom Viewer.

## 9. Personal Conclusion

There is a strong suspicion of malicious behavior in the temporary file associated with the MicroDicom Viewer software, which is not typical for a legitimate image viewer. Further analysis by the software manufacturer is recommended.

## 10. Disclaimer

This report is based on personal analysis and available scan data as of June 2025. It does not constitute a definitive malware diagnosis. Users are advised to perform their own security checks and exercise caution when using the affected software.

---

*Report created by Darius Nasser-Jafar, Cybersecurity Enthusiast*